ekuore

Pulse oximeter

User Manual

EK3011

ekuore

Pulse oximeter

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ekuore

Pulse oximeter

The pulse oximeter is used to check blood oxygen saturation and heart rate. This medical device is reusable and adapts to animals of different shapes and sizes.

It can be used in both operations and postoperative during hospitalization to monitor the patient's condition.

The pack includes:

- eKuore Pulse oximeter [-]
- SpO₂ sensor [-]
- 2 clamps [-]
- Case [-]
- User manual [-]

 CAUTION

- The pulse oximeter is intended only as an adjunct in the patient assessment. It should be used when the patient presents clinical signs and symptoms. The device is not intended to be used as a treatment.
- When using the pulse oximeter together with the electrical surgery equipment, the user must pay attention and ensure the safety of the patient.
- **EXPLOSION HAZARD:** Do not use the pulse oximeter in the presence of flammable anesthetics, explosive substances, vapors, or liquids.
- Make sure not to use the pulse oximeter during MRI (magnetic resonance imaging) scanning or CT (Computed Tomography) environment because induced current could potentially cause burns.
- No modification of the pulse oximeter is allowed. Maintenance should be operated by a professional approved by the manufacturer.
- Please turn off the device before clean the pulse oximeter. Never allow high-pressure and high-temperature disinfection of the device. Never use cleaning agents/disinfectants other than the recommended.
- The pulse oximeter is a sealed product. Keep its surface dry and clean, and prevent any liquid from infiltrating it.
- The pulse oximeter is accurate and fragile. Avoid pressure, knocking, strong vibration, or other mechanical damage. Hold it carefully and gently. When it is not in use, it should be placed in a proper location.
- To dispose of the pulse oximeter and its accessories, follow local regulations or your hospital's policy regarding the disposal of such devices.
- Use AA alkaline batteries. Do not use carbon or poor quality batteries. Remove the batteries if the product will not be used for a long time.
- Please use the SpO2 sensor that corresponds to the product or use one that has been approved by the manufacturer. The user is responsible for checking the compatibility of the pulse oximeter, the sensor, and the cable before use. Incompatible components may cause a degradation in performance. If there are signs of damage, please stop use.
- The user should read the instruction manual carefully and ensure that it is understood or consult the manufacturer before use.
- Avoid static electricity. Before using the pulse oximeter, check the direct or indirect static electricity of all users and patients who will come into contact with the device.
- If the pulse oximeter uses an unspecified and without EMC test system configuration, it can enhance electromagnetic radiation or reduce anti-electromagnetic interference performance. Please use the specified configuration.
- Portable and mobile radio frequency communication equipment can affect the normal use of the pulse oximeter.
- The pulse oximeter should not be near or stacked with other equipment. If this situation occurs during use, check that it works normally with the configuration you are applying.
- Never immerse or soak the pulse oximeter.
- Exercise caution during cleaning/disinfection to avoid wetting the socket.

- Never use cleaning agents/disinfectants other than the recommended.
- Please shut off the power before clean the pulse oximeter. Never permit high-pressure and high-temperature disinfection of the pulse oximeter.

i NOTICE



- Do not place the SpO₂ sensor in the extremities with an arterial catheter or a venous syringe.
- The cable must be connected carefully and prevent the patient from coiling or suffocating.
- Cover the sensor with opaque material in the presence of strong light, otherwise, the measurement will be inaccurate.
- Try to keep the patient still and avoid excessive movement of the measurement point.
- Do not use the pulse oximeter to measure patients whose pulse rate is lower than 30bpm, as this may cause incorrect results.
- Make sure there is no contamination or scarring in the area where the SpO₂ sensor is placed. Otherwise, the results may be incorrect because the signal received by the sensor is affected.
- Disinfection of the pulse oximeter between patients is recommended to avoid cross-contamination.
- Change the position of the sensor after prolonged use or if the patient needs it. After changing the sensor, check the integrity of the skin, the circulatory status, and the correct alignment. It is recommended to make the change every 2 hours.
- Do not perform SpO₂ monitoring and NIBP measurements on the same position simultaneously. Obstruction of blood flow during NIBP measurements may adversely affect the reading of the SpO₂ value.
- The SpO₂ sensor is not allowed to contact the patient's skin if it is above 41°C.
- Shock, anemia, hypothermia, and the application of vasoconstriction drugs may decrease arterial blood flow to an unmeasurable level.
- Pigment or deep color may cause inaccurate measurements.

Device and accessories

1. Up button 2. Down button 3. Function (FN) button 4. Power button*


*Power off: hold the button for 3 seconds

Preparing the device

Assembly

1. Insert the plug of the sensor into the socket of the pulse oximeter.
2. Put the batteries. (Read the "Battery installation" section)

Application to the patient

Install the sensor on the appropriate part of the patient's body (tongue or ear).

Download the App

Download the “eKuore Vet” App for free on your smartphone or tablet.


To download, search for **eKuore Vet App**.

Bluetooth connection

The Bluetooth icon displayed on the screen shows the status of the connection. If the Bluetooth icon:

-  **Not displayed** Bluetooth is enabled.
-  **Flickering** There is no connection to the smartphone/tablet.
-  **Displayed** The pulse oximeter is connected successfully.

First connection

1. Open the "eKuore Vet" app.
2. Press the ON/OFF button to turn it on.

If you have not connected it yet, connect the sensor to the pulse oximeter and the clamp to the patient's body as explained in the "Preparing the device" section.

Make sure that Bluetooth is activated on both devices to get a successful connection.


After opening the app, follow the steps below:

3. Accept the terms of use and privacy policy.
4. Log in with your Apple (iOS) or Google (Android) account.

In the app's main menu, a list of eKuore veterinary devices will be displayed.

5. Find pulse oximeter in the device list and tap on 'Connect'.
6. Accept location permission.
7. Wait until the app finds your device and select pulse oximeter.
8. Now, you can use your pulse oximeter for operations and hospitalizations.

After 10 seconds of inactivity, the device will automatically turn off.

FAQs:

To learn about all the functions of the pulse oximeter, visit:

www.ekuore.com/support

Main functions

Hospitalization


This mode is used to monitor the patient during the postoperative period. It controls blood oxygen saturation levels and heart rate. When the levels are out of range, the app will warn you with alerts previously configured.

Operation

This mode is used during operations, visual and audible alerts will jump, as in hospitalization mode, when levels are outside the configured range.

Also, this mode allows you to record the values of the patient from time to time. This generates a table at the end of the operation that can be shared.

NOTE:

The following graph will be displayed on the screen when the received signal is unsuitable.

Settings

From the device


Hold the FN button for system displays Settings interface. You can set a total of seven parameters.

To browse through the different parameters, press the “FN” button.

Press the “UP” and “DOWN” button to choose the values for each parameter.

1. SpO₂ upper limit* Default value is 99%.
2. SpO₂ lower limit* Default value is 85%.
3. Pulse rate upper limit* Default value is 120bpm.
4. Pulse rate lower limit* Default value is 50bpm.
5. Alert volume
6. Automatic screen out ON/OFF. The screen turns off automatically after 30s of inactivity.
7. Bluetooth ON/OFF.

*When the measured value exceeds the threshold value, the system will remind, the buzzer will be on and the numerical value will flicker.

From the app

You can configure the device from the operation and the hospitalization screen. Press the button on the upper right corner to open the configuration options:

Sound	ON/OFF
Config alarms	SpO ₂ %, tachycardia, bradycardia
Screen	ON/OFF

Pause the alerts

Press the “FN” button for a short time to pause the alerts of the pulse oximeter. The sound will stop but the screen flicker will remain. To set the alert again, just press the “FN” button.

Battery installation

1. Push the locking bar to the right, and then open the battery cover.
2. Install batteries into the slots per the "+" and "-" symbol.
3. Put the battery cover.

You must correctly install the positive and negative electrodes of batteries. Otherwise, the device may be damaged.

When install or remove batteries, please follow the correct steps. Otherwise, the battery compartment may be damaged.


Batteries are not included in the pack.

ELECTRICAL SPECIFICATIONS

Working voltage	D.C.2.32V~D.C.3.4V
Battery type	Two common 1.5V AA alkaline batteries
Power consumption	Smaller than 80mA

Batteries may not be disposed of in household waste (European Battery Law). Please dispose of properly per local ordinance.

Cleaning and disinfection

Cleaning

1. Clean the pulse oximeter with cotton or soft cloth moistened with water.
2. Dry the excess of water with a soft cloth.
3. Allow the pulse oximeter to air dry.

Disinfection

The recommended disinfectants include: 70% ethanol, 70% isopropanol, glutaraldehyde (2%) solution disinfectants.

1. Clean the pulse oximeter as instructed above.
2. Disinfect the pulse oximeter with cotton or soft cloth moistened with one of the recommended disinfectants.
3. After disinfection, be sure to wipe off the disinfectant left on the pulse oximeter with a soft cloth moistened with water.
4. Allow the pulse oximeter to air dry.

NOTE

It is a best practice to use a dry wipe to remove organic material.

Do not use any solvents or IPA, because they can melt the painting or plastic case.

Troubleshooting

TROUBLE	POSSIBLE REASON	SOLUTION
The SpO2 and PR can't be displayed normally and the value disappeared	<ol style="list-style-type: none"> 1. Improper SpO2 sensor installation or incorrect sensor placement on the animal. 2. The animal's SpO2 is too low to be detected. 3. The sensor falls off. 4. Bluetooth signal is interrupted. 	<ol style="list-style-type: none"> 1. Please attach the proper sensor to the proper site on the animal and try again. 2. Try again. 3. Please insert the plug of sensor tight. 4. Check the Bluetooth connection and reconnect.
The SpO2 and PR display instable	<ol style="list-style-type: none"> 1. Sensor placement on the animal is instability. 2. The measured site is shaking or the animal is moving. 	<ol style="list-style-type: none"> 1. Attach the proper sensor to the proper site on the animal and read the SpO2 value after the installation is stability. 2. Let the animal keep calm.
The device can't be powered on	<ol style="list-style-type: none"> 1. The batteries are drained or almost drained. 2. The batteries are not inserted properly. 3. The device's malfunction. 	<ol style="list-style-type: none"> 1. Change batteries. 2. Reinstall batteries. 3. Please contact the supplier.
The screen is suddenly off	<ol style="list-style-type: none"> 1. The product is automatically powered off when no signal is detected longer than 10 seconds. 2. Power quantity of the batteries is exhausted. 	<ol style="list-style-type: none"> 1. Normal. 2. Replace batteries.

Technical features and requirements

APP

OS	Android e iOS
Supported versions	Visit Google Play or App Store to update this information
Requirement	Bluetooth connection enabled

DEVICE

Pulse oximeter size	150mm × 90mm × 26mm
Changer cable length	–
Weight	116g±1g (Do not contain batteries and sensor)
Color	White
Wireless communication	Bluetooth Class II (Distance less than 10m)

More information about warranty at www.ekuore.com/en/warranty
EXCLUSIVE USE FOR VETERINARY FIELD

SpO₂

Measurement range	35~100%
Accuracy	±2% (80%~100%); ±3%(70%~79%)

The method of confirming the blood oxygen measurement accuracy is to compare the oximetry measurement value with the value of blood gas analyzer.

PULSE RATE

Measurement range	25~250bpm
Accuracy	±2bpm

Pulse Rate accuracy has passed proving and comparison with SpO₂ simulator.

ENVIRONMENT REQUIREMENTS

Temperature

Operation	+5~+40°C
Transportation and storage	-10~+50°C

Humidity

Operation	15%~80% (noncondensing)
Transportation and storage	10%~90% (noncondensing)

Atmospheric pressure

Operation	860hPa~1060hPa
Transportation and storage	700hPa~1060hPa

ekuore

Pulse oximeter

ANDROID

iOS

www.ekuore.com

info@ekuore.com

+34 961 156 738

IFU 01_2020/08