# ekuore ECG

Veterinary Monitor 6 Leads

User Manual

# ekuore ECG

## Veterinary Monitor 6 Leads

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# ekuore ECG


## Veterinary Monitor 6 Leads

Digital veterinary electrocardiograph makes it possible to monitor and record patients' cardiac activity wirelessly. Besides, it offers the possibility of digitalising, saving and sending electrocardiograms in PDF to specialists.

### The pack includes:

- eKuore ECG [EK3008]
- 10 disposable electrodes [EK3004-E]
- 3 clamps [EK3004-P]
- 3 Jumpling wires [EK3008-J]
- Quickguide [EK3008-Q]
- Charger cable [EP1500]
- Charger cradle [EK3008-C]
- Case [EK3008-E]

## Device and charging platform

### Status LED

The device's LED light colour code is as follows:

- |                  |                      |
|------------------|----------------------|
| ● Steady red     | Battery charging     |
| ● Flashing green | Device on standby*   |
| ● Steady green   | Monitoring/recording |


\*It appears when returning to the ECG menu after monitoring/recording

## Preparing the device


### Assembly

1. Stick the corresponding stickers that will guide you to take and ECG, as indicated below. Adhesive patches will be applied depending on the standard selected.

#### European standard (IEC-Code)

#### American standard (AHA-Code)

2. Connect the cables to the 6-lead ECG using the clip system included.
3. Insert the male heads into each of the 6-lead ECG female heads.
4. Depending on the patient needs, either patches or clips can be used in the three female head.


**Patches.** The patient must be shaved.

**Clamps.** To use clamps successfully, a conductive liquid such as alcohol or conductive gel must be applied beforehand.


5. Place the electrodes or clips on the patient, following the adhesive patches previously placed as shown in the image.


Patches

Clamp

Cable - Male head

Cable - Female head


**TIP:**

In case of need to take off the clamps, to avoid damaging the device, do NOT PULL the clamps, PUSH SIDEWAYS the clamps


## Application to the patient

European standard

American standard

Download the App

Download the "eKuore Vet" App for free on your smartphone or tablet.

To download, search for eKuore Vet App.

## Device connection

1. Open the "eKuore Vet" app.
2. Hold the ON/OFF button pressed for 3 seconds until the device is turned on (until the green light flashes).


If you have not connected it yet, connect the cables to the 6-lead ECG and the clips or electrodes as explained in the previous section (Preparing the device).

After opening the app, follow the steps below:

3. Accept the terms of use and privacy policy.
4. Log in with your Apple (iOS) or Google (Android) account.

In the app's main menu, a list of eKuore veterinary devices will be displayed.

5. Find 6-lead ECG in the device list and tap on 'Connect'.
6. Accept location permission.
7. Wait until the app finds your device and select 6-lead ECG.
8. After correctly placing the clips or patches, you may now record or monitor using your 6-lead ECG.

### FAQs:

To learn about all the functions of the 6-lead ECG, visit:  
[www.ekuore.com/support](http://www.ekuore.com/support)

## Main functions


### Monitoring

It monitors the patient for an indefinite period of time. It enables sound activation and alarm setting, determining the values for tachycardia and bradycardia. Records are taken in several ranges.

Normal ECG signal

Noise in ECG signal

### Recording

It records the patient during a set period of time. Recording may be stopped at any time. Records are taken in several ranges.

Once the recording ends, the patient's data may be added. All recordings are stored and available for future reference in the App section called Library.

Recordings may be exported in PDF to be shared as a file.


#### TIP:

If the device is correctly connected and placed on the patient, but a loud noise is detected in the ECG signal, other devices in the room may be causing electronic interferences.

## Charge the device

Charge the device connecting the USB cable to the charging platform and let the charge module on the platform.

The charge LED will turn red and when the charge has finished it will turn off.

### BATTERY FEATURES

Rechargeable Lithium Polymer Battery (DC 3.7V, 1000mAh)

Continuous Operating Time 48 hours

Charging Time 2 hours

Power supply 5V 0.5A

# Technical features and requirements

## APP

OS	Android e iOS
Supported versions	Visit Google Play or App Store to update this information
Requirement	Bluetooth connection enabled

## DEVICE

ECG 6-lead size	145mm x 125mm x 12mm
Changer cable length	1m
Weight	47g
Color	White
Wireless communication	Bluetooth Class II (Distance less than 10m)
Consumption	70mW (Transmission signal)

## CLEANING AND STORAGE

It is a best practice to use a dry wipe to remove organic material. Do not use any solvents or IPA, because they can melt the painting or plastic case.

More information about warranty at  
[www.ekuore.com/en/warranty](http://www.ekuore.com/en/warranty)  
**EXCLUSIVE USE FOR VETERINARY FIELD**

Batteries may not be disposed of in household waste (European Battery Law). Please dispose of properly per local ordinance

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